

Kanam

For Barrier Protection and Safety

Providing Solutions to all Your Medical Glove Requirements

surgicare[®]

Low PRO

Your Protection is our Priority

Latex Prepowdered Surgical Gloves with protein content less than $100\mu\text{g}/\text{dm}^2$

- ▶ Latex protein content less than $100\mu\text{g}/\text{dm}^2$ to reduce risk of user being sensitized to Latex Protein Allergy Type I.
On special request Protein Content less than $50\mu\text{g}/\text{dm}^2$ can be provided.
- ▶ Non-detectable levels of chemical residue reduces allergic skin reactions like Contact Dermatitis - Type IV Allergy
- ▶ Lightly powdered with Modified Bio-Absorbable Corn starch USP Grade, to prevent wound infections and granulomas.
- ▶ Manufactured in clean environment to achieve bio burden level of less than 20 cfu prior to Sterilization.
- ▶ Sterilized in a validated process cycle to achieve sterility assurance level of 10^6 , and air wash cycles to reduce sterilitant residue to less than 5 ppm.
- ▶ Latex Film structure passes Viral Penetration Test (ASTM F1671.976)
- ▶ AQL < 1.5 for pin holes
- ▶ Passes Bio-Compatibility Test as per ISO 10993-10

An ISO 9001-2000, ISO 13485 & ISO 14001 Certified Company Inspected by, and Registered with USFDA 510K & Bureau of Indian Standards ISI 13422 & ISO 10282

Micro rough surface for better instrument grip without compromising on sensitivity

Thumb ball design to ease stress on thumb during surgical procedures

Low Allergenicity

Meets standards of
ASTM D-3577
EN 455-1,2 & 3
ISI 13422 &
ISI 4148

Precision designed fit for optimum comfort

CE
0434

SURGICARE
Low Pro



Kanam Latex Industries Pvt. Ltd.
Nagercoil, India
Fax: 91-481-2563614
e-mail: kanam@asianetindia.com
www.kanamlatex.com