

# Kanam

Providing Solutions to all Your Medical Glove Requirements

## surgicare®

Your Protection is our Priority

**Latex Prepowdered Surgical Gloves**

- ▶ Non-detectable levels of chemical residue reduces allergic skin reactions like Contact Dermatitis - Type IV Allergy
- ▶ Lightly powdered with Modified Bio-Absorbable Corn starch USP Grade (No talc), to prevent wound infections and granulomas.
- ▶ Manufactured in clean environment to comply with EN 556 and ISO 13485 requirements.
- ▶ Sterilized in a validated process cycle to achieve sterility assurance level of  $10^6$ , and air wash cycles to reduce sterilant residue to less than 5 ppm.
- ▶ Latex Film structure passes Viral Penetration Test (ASTM F1671.976)
- ▶ AQL 1.5 for pin holes

Precision designed fit for optimum comfort


**CE 0434**

An ISO 9001-2000, ISO 13485 & ISO 14001 Certified Company Inspected by, and Registered with USFDA 510K & Bureau of Indian Standards ISI 13422 & ISI 4148

*It's not just another glove!*

Micro rough surface for better instrument grip without compromising on sensitivity

Meets standards of  
ASTM D-3577  
EN 455-1,2 & 3  
ISI 13422 &  
ISI 4148



Manufactured within the quality guidelines of ISO 13485 and EN 455



**Kanam Latex Industries Pvt. Ltd.**  
Nagercoil, India  
Fax: 91-481-2563614  
e-mail: [kanam@asianetindia.com](mailto:kanam@asianetindia.com)  
[www.kanamlatex.com](http://www.kanamlatex.com)

SURGICARE