



Powderfree Sterile & Non Sterile Latex Examination Gloves

Kaltex Powderfree Sterile and Non Sterile Latex Examination Gloves are made from Natural Rubber Latex. These gloves provide excellent barrier protection for both healthcare worker and patient.

SAFETY FEATURES

- ▶ Enhanced elasticity, tear strength and donning of powderfree gloves.
- ▶ Beaded cuff with minimum length of 240 mm. Can offer 290 mm length also.
- ▶ Sterilized in a validated process cycle as per ISO 11135:2014 norms for sterile products for sterile examination gloves.
- ▶ Latex protein no more than 50µg of protein per dm², eliminates risk of Latex Allergy - Type 1.
- ▶ Absence of powder eliminates dangers associated with powder, preventing granulomas, adhesions, fibrosis and occupational asthma. Does not contaminate sterile work area.

STANDARDS

- ▶ EN 455 -1: 2000, EN 455 -2: 2015, EN 455 -3; 2015 & ASTM D-3578: 2009.
- ▶ CE Mark - MDD Class Is/PPE Cat. III.
- ▶ Registered with USFDA (510K).

QUALITY CONTROL

- ▶ Factory standard of 1.5 AQL for pinholes.



Kanam
KANAM LATEX
INDUSTRIES PVT. LTD.

The trusted name the world over

Ooppottil Buildings, K.K. Road,
Kottayam-686 001, Kerala, India
Tel: 0481-2560108, 2563513
Fax: 0481-2563614
E-mail: customer@kanamlatex.com
kanamlatex@gmail.com

www.kanamlatex.com

POWDER FREE STERILE & NON STERILE LATEX EXAMINATION GLOVES

Product Information

TYPE

Powderfree Sterile & Non Sterile Latex Examination Gloves

PRIMARY MATERIAL

Natural Rubber Latex (NR).

COLOUR

White to Pale Yellow.

DESIGN

Ambidextrous, Straight Finger, Beaded Cuff,
Micro Rough Textured.

POWDER FREE RESIDUE

<2mg/glove.

STERILIZATION

Ethylene Oxide for sterile examination gloves.

SHELF LIFE

3 years in original packaging if stored as per storage instructions.

PHYSICAL PROPERTIES

Characteristics	Before Ageing	After Ageing 70±2°C for 166±2hrs
Tensile Strength (Mpa)	22 min	18 min
Ultimate Elongation (%)	800 min	600 min
Force at break (N) Min.	6 During Shelf Life (N)	6 After Challenging (N)

PHYSICAL DIMENSIONS

Size	Length (mm) Min	Palm Width (mm)	Average weight per pair (gm)	Thickness (mm) Min	
				Palm	Finger
Small	240	84 - 85	5.5 - 6.0	0.10	0.11
Medium	240	95 - 96	6.0 - 6.5		
Large	240	106 - 107	6.5 - 7.0		
X-Large	240	112 - 113	7.0 - 7.5		

PACKAGING (in cm) - Sterile


Material	Quantity	Dimension
Paper Wallet(Folded)	One Pair	15 (L) x 9.5 (W)
Pouch	One Pair	20 (L) x 13 (W)
Inner Box	50 Pairs	20(L) x 13(W) x 18(H)
Master Carton	300 Pairs	40(L) x 21(W) x 37(H)

PACKAGING (in cm) - Non Sterile

Material	Quantity	Dimension
Inner Box	100 Pcs	23.8(L)x13(W)x7.2(H)
Master Carton	1000 Pcs	38(L)x27(W)x25(H)



Powder Free Sterile & Non Sterile Latex Examination Gloves



KLIP/IG/17/Rev#02