

# EU TYPE-EXAMINATION CERTIFICATE


## Regulation (EU) 2016/425, MODULE B

### 0598/PPE/25/3450 Issue 1

<b>Product</b>	Protective gloves against dangerous chemicals and micro-organisms
<b>Model</b>	Powder Free Dual Latex Surgical Gloves With Puncture Indication (Sterile & Non-Sterile)
<b>Trademark</b>	Surgicare Dual, Dermagel Dual, Surgicare Dual Premium
<b>Certificate Holder / Manufacturer</b>	Kanam Latex Industries Private Limited 12/67 C, Ananthanadarkudy, Asaripallam P.O, Nagercoil – 629 201, Kanyakumari District, Tamil Nadu, India
	Product complies with the applicable essential health and safety requirements of Regulation (EU) 2016/425 and standard(s) mentioned below
<b>Standard(s)</b>	EN ISO 374-1:2016+A1:2018 (Protective gloves against dangerous chemicals and micro-organisms) Type B glove against reagents: Under glove: K, M, P, T, Over glove: K, M, T; EN ISO 374-5:2016 (Protective gloves against dangerous chemicals and micro-organisms) gloves protecting against bacteria and fungi (EN 374-2:2019 inspection level 2 (AQL <1.5 G1), and virus (ISO 16604)); EN ISO 21420:2020 (Protective gloves - General Requirements).
<b>Other Information</b>	This certificate shall be used in conjunction with conformity assessment procedure module C2 or D.
<b>Validity</b>	This certificate is valid until 2030-05-27.
<b>Date of issue</b>	2025-05-27

**SGS Fimko Oy**

**Signature**



Kelvin Shepherd

SGS Fimko Oy is a Notified Body (0598) according to the Personal Protective Equipment Regulation (EU)



**Additional information** Powder Free Dual Latex Surgical Gloves With Puncture Indication (Sterile & Non-Sterile), two different products in a single pack:

Overglove: Powder Free Polymer Coated Latex Surgical Gloves (Sterile & Non-Sterile), white to pale yellow natural rubber latex glove, anatomic shape, beaded cuff, polymer coated, single use, available in eight nominal hand sizes: 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5 & 9.0.

Underglove: Powder Free Latex Surgical Gloves as an Underglove (Sterile & Non-Sterile), green natural rubber latex glove, anatomic shape, beaded cuff, polymer coated, single use, available in eight nominal hand sizes: 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5 & 9.0.

Sold under brand names: Surgicare Dual, Dermagel Dual, Surgicare Dual Premium.

**Contract No:** INP-CERT250100006-01/GB-202410007406.

**Technical File:** TCF/PPE/DLG-PF, Rev 00, Dated: 05.05.2025.