

# EU TYPE-EXAMINATION CERTIFICATE

## Regulation (EU) 2016/425, MODULE B

### 0598/PPE/25/2382 Issue 1

**Product** Protective gloves against dangerous chemicals and micro-organisms  
**Model** Powder Free Latex Surgical Gloves (Sterile & Non-Sterile)

**Trademark** Surgicare Plus, Comfort PF, Comfort Powder-Free

**Certificate Holder / Manufacturer** Kanam Latex Industries Private Limited  
 12/67 C, Ananthanadarkudy, Asaripallam P.O, Nagercoil – 629 201, Kanyakumari District, Tamil Nadu, India

Product complies with the applicable essential health and safety requirements of Regulation (EU) 2016/425 and standard(s) mentioned below

**Standard(s)** EN ISO 374-1:2016+A1:2018 (Protective gloves against dangerous chemicals and micro-organisms) Type B glove against reagents: K, M, P & T;  
 EN ISO 374-5:2016 (Protective gloves against dangerous chemicals and micro-organisms) gloves protecting against bacteria and fungi (EN 374-2:2019 inspection level 2 (AQL <1.5 G1), and virus (ISO 16604));  
 EN ISO 21420:2020 (Protective gloves - General Requirements).


**Other Information** This certificate shall be used in conjunction with conformity assessment procedure module C2 or D.

**Validity** This certificate is valid until 2030-02-12.

**Date of issue** 2025-02-12


**SGS Fimko Oy**

**Signature**



Kelvin Shepherd

SGS Fimko Ltd is a Notified Body (0598) according to the Personal Protective Equipment Regulation (EU)



**Additional information** Powder Free Latex Surgical Gloves (sterile & non-sterile), white to pale yellow natural rubber latex glove, anatomically shaped, beaded cuff, textured external surface, chlorinated internal surface, single use, which are available in eight nominal hand sizes: 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5 & 9.0.

Sold under brand names: Surgicare Plus, Comfort PF, Comfort Powder-Free.

**Contract No:** INP-CERT250100002-01/GB-202410007359.

**Technical File:** TCF/PPE/LSG-PF, Issue 02, Rev 06, Dated: 23.10.2024.